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MAJOR CAR SEAT

INSTRUCTION MANUAL

This child restraint complies with ECE regulations for use by children:

- Group I / II / III
- Mass: 9 – 18 kg with built-in harness
- Mass: 15 – 36 kg as high-back booster

READ THIS MANUAL

Do not install or use this child restraint until you read and understand the instructions in this manual and in your vehicle owner's manual.

FAILURE TO PROPERLY USE THIS CHILD RESTRAINT INCREASES THE RISK OF SERIOUS INJURY OR DEATH IN A SHARP TURN, SUDDEN STOP OR CRASH.


IF YOU HAVE ANY PROBLEMS, PLEASE DO NOT RETURN THIS PRODUCT TO THE PLACE OF PURCHASE. PLEASE CONTACT CUSTOMER SERVICE AT +49 4761 886 0 OR AT WWW.THOMASHILFEN.COM

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If you need help.....	3
Warnings	3
General warnings	3
Usage warnings	5
Injury warnings	6
Vehicle warnings	7
Certifications	7
Airbag warnings	8
Damage warnings	9
Mass and height limitations	10
With built-in harness	10
Without built-in harness	11
Maximum height	11
Parts list	12
Vehicle seating positions	14
Safe vehicle seat belt systems	15
Unsafe vehicle seat belt systems	16
Assembling/disassembling backrest	17
Preparing child restraint for use	18
Adjusting harness shoulder straps	18
Headrest height adjustment.....	20
Installing with vehicle seat belt	20
Using child restraint with built-in harness	22
Securing child in child restraint	22
Removing child from child restraint	24
Using child restraint in booster mode (without internal harness)	24
Removing/installing harness	24
Removing the crotch buckle	26
Preparing booster seat for use	27
Securing child and booster seat in vehicle	28
Installation - booster seat	29
Using cup holder (when applicable)	31
Aircraft installation	31
Care and maintenance	31

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IF YOU NEED HELP

Please fill in the information below for future reference. The model number and manufacturing date can be found on a label on the underside of the child restraint.

Model Number: _____

Manufacture Date: _____

Date of Purchase: _____

Our Customer Service Department is standing by to help you with any questions or concerns you may have regarding installation, use, care and/or parts. Please have your product's model number and manufacturing date available when contacting us.

WARNINGS

GENERAL WARNINGS

WARNING!

The safety of your child and the effectiveness of this child restraint depend on:

- 1. Selecting a suitable location** for the child restraint in your vehicle.
- 2. Properly routing the vehicle seat belt.**
- 3. Properly securing both your child and the child restraint in the vehicle** using a seat belt designed to restrain both the child restraint and your child. Many seat belts are **NOT** safe to use with this child restraint. Please review the sections "Safe/Unsafe vehicle seat belt systems" carefully. If in doubt about installing this child restraint, consult the vehicle owner's manual or call the Thomashilfen Customer Service Department.

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⚠ WARNING!

FAILURE TO PROPERLY READ, FOLLOW AND USE THESE INSTRUCTIONS AND THE CHILD RESTRAINT LABELS INCREASES THE RISK OF SERIOUS INJURY OR DEATH IN A SHARP TURN, SUDDEN STOP OR CRASH. Your child's safety depends on the proper installation and use of this child restraint. Even if this child restraint seems easy to install and use, it is **VERY IMPORTANT** to read and understand these instructions and the vehicle owner's manual.

Follow all labels on this child restraint and in these written instructions. **DO NOT INSTALL OR USE THIS CHILD RESTRAINT UNTIL YOU READ AND UNDERSTAND THE INSTRUCTIONS IN THIS MANUAL AND IN YOUR VEHICLE OWNER'S MANUAL.**

NO CHILD RESTRAINT CAN GUARANTEE PROTECTION WHEN THE VEHICLE IS SERIOUSLY IMPACTED. However, when correctly installed, a child restraint will substantially improve the chances for survival and reduce the risk of serious injury in most crashes.

Carefully follow the instructions in this manual and in your vehicle owner's manual. If there is a discrepancy between the two, the vehicle owner's manual must take precedence with regards to seat belts and/or airbags.

! ADULT ASSEMBLY REQUIRED.

! USE THIS CHILD RESTRAINT FORWARD-FACING ONLY.

! DO NOT USE THIS CHILD RESTRAINT without properly assembling it.

! DO NOT allow children to play with this child restraint.

! NEVER leave children alone in a vehicle under any circumstance.

! It is recommended that children are not left in their child restraint system unattended.

! Set a good example by ALWAYS wearing your seat

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belt. Make it a firm rule that the vehicle does not move until everyone is buckled up. **MAKE NO EXCEPTIONS.**

! It is not advised to loan or pass along a child restraint. **DO NOT USE A SECOND-HAND CHILD RESTRAINT** or a child restraint whose history you do not know.

! When disposing a child restraint system, ensure to dispose of the unit in a way others cannot reuse it (e.g. use a dark trash bag, do not sell to a thrift shop or in a garage sale, etc.).

! Store the child restraint in a safe place when not in use. **DO NOT** place heavy objects on top of the child restraint while in storage.

! Do not remove any stickers or warnings as they contain important information.

USAGE WARNINGS

WARNING!

! **ALWAYS** ASSIST YOUR CHILD WHEN GETTING INTO AND OUT OF THE CHILD RESTRAINT.

! **NEVER** LET YOUR CHILD LEAN ON THE ARMREST WHEN GETTING INTO OR EXITING THE BOOSTER SEAT TO PREVENT THE SEAT FROM TIPPING.

! **ENSURE** THAT THE CHILD RESTRAINT IS LEVEL, STABLE AND DOES NOT ROCK SIDE-TO-SIDE. IF THE CHILD RESTRAINT IS NOT STABLE, TRY ANOTHER SEATING LOCATION.

! **DO NOT ADD ANYTHING BEHIND OR UNDER THE CHILD RESTRAINT.**

According to accident statistics, **CHILDREN ARE SAFER WHEN PROPERLY RESTRAINED IN THE REAR SEATING POSITIONS THAN IN THE FRONT SEATING POSITIONS.**

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Please study the section "Vehicle seating positions".

! **DO NOT** use any load bearing contact points other than those described in the instructions and marked in the child restraint.

! **DO NOT** allow others to install or use this child restraint unless they read and understand all the instructions and labels.

! **DO NOT GIVE THIS CHILD RESTRAINT** to someone else without also giving them this instruction manual. Additional instruction manuals can be downloaded directly from www.thomashilfen.com.

! To ensure proper child restraint use, keep and store this instruction manual and refer to it when necessary.

! **DO NOT LEAVE A CHILD RESTRAINT UNSECURED IN YOUR VEHICLE.** In a collision, sharp turn or sudden stop, an unsecured child restraint can injure other occupants. Remove it or make sure that it is securely attached with the vehicle seat belt, even when not in use.

INJURY WARNINGS

! **INJURIES MAY OCCUR EVEN IF YOU ARE NOT IN A CRASH** if the child restraint is not properly installed or if your child is not properly secured in the child restraint.

! The rigid items and plastic parts of a child restraint shall be so located and installed that they are not liable, during everyday use of the vehicle, to become trapped by a movable seat or in a door of the vehicle.

! **DO NOT** leave folding vehicle seats unlatched. In the event of a sharp turn, sudden stop or crash, a loose seat back could prevent the child restraint from protecting the child as well as it should.

! It is recommended that any luggage or other objects liable to cause injuries in the event of a collision shall be properly secured. In the event of a sharp turn,

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sudden stop or crash, loose objects will keep moving, potentially causing serious injuries.

! The child restraint shall not be used without the cover. It **MUST** always be used with the fabric cover provided by the manufacturer.

! The seat cover should not be replaced with any other than the one recommended by the manufacturer, because the cover constitutes an integral part of the restraint performance.

! **PARTS OF THE CHILD RESTRAINT COULD BECOME HOT ENOUGH TO BURN THE CHILD IF LEFT IN THE SUN.** Contact with these parts can injure your child. **ALWAYS** touch the surface of any metal or plastic parts before putting your child in the child restraint. Cover the child restraint with a towel or blanket when the vehicle is parked in direct sunlight. **KEEP AWAY FROM SUNLIGHT, OTHERWISE IT MAY BE TOO HOT FOR THE CHILD'S SKIN.**

VEHICLE WARNINGS

! **DO NOT** use this child restraint with 2-point vehicle belts (lap-belt only or shoulder-belt only). Your child would be seriously injured in the event of a crash.

! **DO NOT** use this child restraint with unsafe vehicle seat belts listed in section "Unsafe vehicle seat belt systems" of this manual.

! **DO NOT** use this child restraint if it can't be properly installed. It is possible that your vehicle has no suitable seating position for this child restraint due to seat belt type, airbag system and/or lack of location of a tether strap anchor. Read and follow these instructions, and your vehicle owner's manual carefully.

CERTIFICATIONS

This is a "Universal" child restraint when installed using the 3-point vehicle seat belt only. It is approved to the

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requirements of the European Standard for Child Safety Seat Equipment under Regulation No. 44, 04 series of amendments, for general use in vehicles and it will fit most, but not all, car seats.

A correct fit is likely if the vehicle manufacturer has declared in the vehicle handbook that the vehicle is capable of accepting a "Universal" child restraint for this age group.

This child restraint has been classified as "Universal" under more stringent conditions than those which applied to earlier designs which do not carry this notice.

If in doubt, consult either the child restraint manufacturer or the retailer.

! Only suitable for use in the listed vehicles fitted with 3-point with retractor safety-belts, approved to ECE Regulation No. 16 or other equivalent standards.

This restraint is **NOT CERTIFIED** for aircraft use.

This child restraint must **NOT** be used in homes, boats, buses, motorcycles, motor homes, or any other vehicles that have non-forward-facing seats, missing seat belts, or whose seat belts do not meet the requirements of this child restraint.

AIRBAG WARNINGS

! **NEVER PLACE THIS CHILD RESTRAINT IN A VEHICLE SEATING LOCATION THAT HAS AN ACTIVE FRONTAL AIRBAG.** An inflating airbag may cause serious injury or death. Refer to your vehicle owner's manual for information about airbags and child restraint installation.

This child restraint may be used in a rear seating position equipped with side impact airbags providing they are:

- Roof mounted;
- "Curtain" type deployment;
- Will **NOT** interfere with, impact or in any way reposition the child restraint when inflated.

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! If a vehicle's front seat location must be used, **ALWAYS** disable the airbag. Failure to disable the airbag system can result in serious injury or death. Follow the vehicle owner's manual for instructions on disabling the front airbag system.

Please refer to your vehicle owner's manual for airbag type, and for appropriate use with child restraints.

DAMAGE WARNINGS

! It is recommended that the device should be changed when it has been subject to violent stresses in an accident.

! **DO NOT USE A CHILD RESTRAINT THAT HAS BEEN IN AN ACCIDENT OF ANY KIND.** An accident can cause damage to the child restraint that you cannot see.

! **DO NOT** use this child restraint if it is damaged in any way or has missing parts.

! **NEVER** use a cut, frayed or damaged vehicle seat belt or harness.

! **DO NOT** tamper with or modify the child restraint, for any reason. Any unauthorised or additions to the device without the approval of the Type Approval Authority can cause serious injury or death in a sharp turn, sudden stop, or crash.

! **DO NOT** use anything to raise the child restraint off the vehicle seat. Untested modifications, such as adding extra padding, could affect the safety of the child restraint and result in serious injury or death in a crash.

! **DO NOT, EXCEPT AS DESCRIBED IN THESE INSTRUCTIONS, MODIFY YOUR CHILD RESTRAINT, DISASSEMBLE ANY PART OR CHANGE THE WAY THE VEHICLE SEAT BELTS ARE USED.** This includes any accessories or parts by other manufacturers as non-tested accessories can affect the safety of the child restraint.

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⚠ WARNING!

THIS CHILD RESTRAINT HAS DIFFERENCES IN INSTALLATION AND USAGE DEPENDING ON THE WEIGHT AND HEIGHT OF YOUR CHILD. FAILURE TO USE THIS CHILD RESTRAINT WITHIN THESE GUIDELINES MAY INCREASE THE RISK OF SERIOUS INJURY OR DEATH.

! WARNING: Verify your child's height and weight frequently. **DO NOT** guess before choosing the proper setup for your child.

It is recommended that you **ALWAYS** confirm country laws with local authorities. Laws vary from country to country regarding minimum/maximum child restraint requirements for age, mass and height.


WITH BUILT-IN HARNESS

To use this child restraint, your child **MUST MEET ALL** of the following requirements:

• GROUP I

Forward-facing with built-in harness:
between 9 - 18 kg
(20 - 39.6 lbs.)

- Shoulders are **BELOW** top harness slots ①.
- Child's ears **MUST BE BELOW** top of the headrest. If tops of ears are above the top of the headrest at its maximum height, the child is too large for the child restraint ②.

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To use this child restraint, your child **MUST MEET ALL** of the following requirements:

• **GROUP II / III
Forward-facing
with backrest:**

between 15 - 36 kg
(33 - 79 lbs.)

- Shoulder belt **MUST** lie across the **MIDDLE** of the child's shoulder and chest, and lap portion **MUST** be positioned low across child's hips ①.
- Child's ears **MUST BE BELOW** top of booster seat headrest at its maximum height. If tops of ears are above the top of the headrest at its maximum height, the child is too large for the booster seat ②.

MAXIMUM HEIGHT

⚠ WARNING!

Child is too tall when top of ears (mid-point of child's head) are above the top of the child restraint/booster seat headrest at its maximum height.

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- A. Headrest
- B. Backrest
- C. Base
- D. Shoulder belt positioning guides
- E. Harness straps, buckle and clips
- F. Shoulder pads
- G. Crotch strap and pad
- H. Metal up-front harness adjustment button
(under fabric seat pad)
- I. Up-front harness adjustment strap
(under fabric seat pad)
- J. Seat belt installation guides
(under fabric, when applicable)
- K. Cup holder (when applicable)
- L. Head insert cushion (when applicable)
- M. Body insert cushion (when applicable)
- N. Headrest height adjustment button
- O. Harness splitter plate
- P. Shoulder pad splitter plate

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

TO PURCHASE REPLACEMENT PARTS:

In Germany, please contact us at
www.thomashilfen.de
 or by calling +49 4761 886 0.

under fabric

⚠ WARNING!

FAILURE TO PROPERLY SECURE YOUR CHILD AND THE CHILD RESTRAINT IN YOUR VEHICLE INCREASES THE RISK OF SERIOUS INJURY OR DEATH IN THE EVENT OF A SHARP TURN, SUDDEN STOP OR CRASH.

⚠ WARNING!

IMPROPER PLACEMENT OF THE CHILD RESTRAINT INCREASES THE RISK OF SERIOUS INJURY OR DEATH.

Some seating positions in your vehicle may **NOT** be safe for this child restraint. Some vehicles have **NO** seating positions which are compatible with this child restraint due to seating surface, seat belt type and/or airbag system. Read and follow these instructions, and your vehicle owner's manual carefully.

⚠ WARNING!

Children aged 13 and under **MUST** ride in the rear seats, away from front airbags. The only exception is if the airbags are turned off. When using a child restraint in a vehicle with active airbags, refer to the vehicle owner's manual for additional installation instructions and precautions. For rear seat installation, adjust the front seat forward enough to ensure the child's feet do not hit the backrest of the front seat.

⚠ WARNING!

Only use this child restraint in forward-facing vehicle seats. Side-facing or rear-facing seats in vans, station wagons or trucks **CANNOT** be used.

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THE VEHICLE SEAT MUST FACE FORWARD. DO NOT use the child restraint with the following vehicle seats:

- Seats which face sides of vehicle.
- Seats which face rear of vehicle.
- Rear seats of small extended cab pickup trucks.

- The front of the child restraint **MUST NOT** hang over the edge of the vehicle seat by more than 5 cm (2 inches).

Refer to the vehicle owner's manual for child restraint installation instructions for your specific vehicle.

! DO NOT use this child restraint in a seating position where the vehicle seat belt or seat belt buckle is not easily accessible with the child restraint secured in place.

SAFE VEHICLE SEAT BELT SYSTEMS

WARNING!

SEAT BELTS DIFFER FROM VEHICLE TO VEHICLE AND FROM SEATING LOCATION TO SEATING LOCATION WITHIN THE SAME VEHICLE. NOT ALL VEHICLE SEAT BELTS CAN BE USED WITH THIS CHILD RESTRAINT.

You must determine if your vehicle's seat belts are compatible for use with this child restraint. If you are unsure, consult your vehicle owner's manual or call Thomashilfen.

! This child restraint is suitable for use in the listed vehicles fitted with 3-point with retractor safety-belts (lap-and-shoulder), approved to ECE Regulation No. 16 or other equivalent standards.

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

Review your vehicle owner's manual for vehicle seat belt use with child restraints.

UNSAFE VEHICLE SEAT BELT SYSTEMS

This child restraint can be installed properly using **MOST** vehicle 3-point lap-and-shoulder seat belts. However, some vehicle seat belts are not safe to use as their anchoring points for the lap belt may not allow for a secure installation. **NEVER** use any unsafe vehicle seat belt systems.


DO NOT USE ANY OF THE FOLLOWING VEHICLE SEAT BELT SYSTEMS WITH THIS CHILD RESTRAINT:

! NEVER use a lap-only belt.


! NEVER use a non-locking vehicle seat belt.

! NEVER use a cut, frayed or damaged vehicle seat belt.

! NEVER use motor driven automatic seat belts.

! NEVER use a vehicle seat belt system where the lap and shoulder belts have independent retractors.

! NEVER use a vehicle seat belt system with top or bottom anchorage points attached to the door in any way.

! NEVER use lap belts whose anchor points are forward of the seat crease.

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


1. Insert 'U-Shape' plastic hooks on bottom of backrest into the corresponding bars on seat base ①.
2. Apply a small amount of force to gently lock the hooks onto the base ②.
3. Slowly swing the backrest upwards, whereby locking the backrest into place ③.
4. To remove backrest, repeat previous steps in reverse order.

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
! WARNING!

INCORRECT USE OF THE HARNESS SYSTEM IS THE MOST COMMON MISTAKE PEOPLE MAKE. IT IS VERY IMPORTANT TO READ, UNDERSTAND AND FOLLOW ALL INSTRUCTIONS AND WARNINGS. DO NOT MISTAKE COMFORT FOR SAFETY. THE HARNESS SYSTEM MUST BE SNUGLY ADJUSTED AT ALL TIMES.

! WARNING: The shoulder harness must be **EVEN OR SLIGHTLY ABOVE** the child's shoulders. If the upper most harness position is below the child's shoulders, **DO NOT** use this child restraint with harness.



! ENSURE the harness height is correct for your child's size. **VERIFY OFTEN.**

! WARNING: DO NOT use this child restraint if the midpoint of your child's head is above the top of the child restraint's headrest at its maximum height.

ADJUSTING HARNESS SHOULDER STRAPS

1. Loosen harness system by pushing the up-front metal harness adjuster at the front of the seat just above the adjustment strap (under the fabric) **1** and by simultaneously pulling the chest strap towards you **2**.

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
TIP: PULL ON THE HARNESS STRAPS. DO NOT PULL ON THE HARNESS PADS AS THEY ARE LOCKED IN PLACE.

2. Face the rear of the child restraint and remove the harness straps and shoulder straps from the splitter plates (see "Removing/installing harness") ③.


3. Gently pull the straps out from the front harness slots.
4. Select the proper harness slot height and feed the straps (harness and shoulder pads) back through. **ENSURE THAT THE HARNESS STRAPS ARE THREADED AT THE SAME LEVEL, ARE NOT TWISTED AND ARE IN THE CORRECT SLOTS.**
5. Reattach the harness straps and shoulder straps to the metal splitter plates.

! ENSURE BOTH HARNESS STRAPS ARE PROPERLY SECURED ON THE SPLITTER PLATE.

WARNING!

FAILURE TO REINSTALL THE HARNESS STRAPS CORRECTLY COULD RESULT IN SERIOUS INJURY OR DEATH.

REMEMBER to check the shoulder straps position frequently and move them to the next position when necessary. The harness straps **MUST** be in the height position **EVEN WITH** or just **ABOVE** the shoulders. The harness pads and crotch pads are optional and are not a safety feature.

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
Your child restraint has seven (7) height adjustments.


- Use one hand to press the headrest adjustment button found on the back of the seat ① while simultaneously pulling the headrest up or down to the desired height with the other hand ②.

Adjust headrest height so that the bottom of the headrest is **EVEN** with the top of the child's shoulders ③.

INSTALLING WITH VEHICLE SEAT BELT

1. Place the child restraint flush against a forward-facing vehicle seat. The vehicle seat **MUST** be equipped with a 3-point lap-and-shoulder belt. The child restraint must be level, stable and not rock side-to-side (see "Vehicle seating positions").
2. Loosen the harness system. Pass the vehicle seat belt behind the harness system over the base of the unit as shown and buckle the vehicle seat belt. Listen for an audible "click". The vehicle seat belt **MUST NOT** be twisted.
3. Lift the back portion of the Base Seat Pad and place the lap portion of the vehicle seat belt **FLAT** against the plastic base ①. **ENSURE**

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

that the shoulder portion of the vehicle seat belt is **FULLY INSERTED IN THE GROOVE BETWEEN THE BACKREST AND BASE** and the lap portion is placed under the armrest ②.

4. While lifting the backrest pad, guide the shoulder belt into the Seat Belt Installation Guide (Parts list - 'J'). See illustrations for proper routing path of the 3-point lap-and-shoulder belt ③.

5. If your vehicle is equipped with a ELR/ALR switchable seat belt retractor system, we recommend that you pull out the full length of the vehicle seat belt (until you cannot pull out any more slack). This locks the retractor. If not equipped with switchable ELR/ALR retractor, continue with step #6.

6. Tighten vehicle seat belt by placing your knee on the child restraint base and use your full weight to push the child restraint down ④ as you simultaneously release the extra seat belt length until seat belt is taught and has no slack ⑤.

For more information on how to tighten the seat belt for your seat belt system, see your vehicle owner's manual.

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WARNING!

FAILURE TO FASTEN AND TIGHTEN THE HARNESS SYSTEM CORRECTLY MAY CAUSE THE CHILD TO BE EJECTED FROM THE RESTRAINT SYSTEM IN CASE OF A SHARP TURN, SUDDEN STOP OR CRASH CAUSING SERIOUS INJURY OR DEATH. In a crash, the child restraint will provide **NO** protection for your child unless the harness straps are properly secured. **DO NOT** mistake comfort for safety. Harness system **MUST** be snugly adjusted. Follow these instructions carefully.

! Incorrect use of the harness system is the most common mistake people make. It is **VERY** important to read, understand and follow **ALL** instructions and warnings. **DO NOT** mistake comfort for safety. The harness system **MUST** be snugly adjusted at all times.

! Snugly adjust the harness provided with this child restraint around your child. A snug strap is one with no slack. It should not press on the child's flesh or cause your child's body to be in an unnatural position. You should only be able to slide **ONE FINGER** between the shoulder straps and the child's shoulder.

! Avoid unnecessary bulky clothing.

SECURING CHILD IN CHILD RESTRAINT

1. To loosen the harness strap: Push the up-front metal harness adjuster at the front of the seat just **ABOVE** the adjustment strap (under the fabric) **1** and simultaneously pull the chest strap towards you **2**.

TIP: PULL ON THE HARNESS STRAPS. DO NOT PULL ON THE HARNESS PADS AS THEY ARE LOCKED IN PLACE.

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2. Unbuckle harness straps by pressing on the **RED** buckle release button and by pulling the buckle clips out **3**. Move harness straps to the side.






3. Place child in child restraint with child's back flat against back of the backrest.

4. Adjust headrest height so that the bottom of the headrest is even with the top of the child's shoulders.

5. Place harness straps over child's shoulders and insert the two-piece buckle clip into buckle until they 'click' into place **4**.

! ENSURE that the buckle is properly fastened by gently tugging on it and that the colour indicator has turned to **GREEN** **5**.

! ENSURE that the straps are **NOT** twisted.

! DO NOT use child restraint if harness does not lock.

6. Tighten harness system by gently pulling on the up-front harness adjuster strap until all slack is removed from the system **6**.

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
! ENSURE straps are flat and snug against child's shoulders and thighs. You should **NOT** be able to slide more than one finger between the strap and the child's shoulder **7**.

A snug strap should not allow any slack. It lies in a relatively straight line without sagging. It should not press the child's flesh or cause your child's body to be in an unnatural position.

! Any straps holding the restraint to the vehicle should be tight, any straps restraining the child should be adjusted to the child's body, and all straps should not be twisted.

REMOVING CHILD FROM CHILD RESTRAINT

1. Loosen the harness system by pressing on the up-front metal harness adjuster (under the fabric), and simultaneously pulling the shoulder straps towards you.
2. Unbuckle harness by pressing on the **RED** buckle release button.

USING CHILD RESTRAINT IN BOOSTER MODE (WITHOUT INTERNAL HARNESS)

REMOVING/INSTALLING HARNESS

1. Keep the harness buckled. Loosen the harness system by pressing on the up-front metal harness adjuster (under the fabric), and simultaneously pulling the shoulder straps towards you.

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
2. Face the rear of the child restraint and remove the harness straps and shoulder straps from the splitter plates ①.


3. Pull the harness straps out of their respective slots in the backrest ② and remove shoulder pads.


4. Continue to thread the harness strap through the plastic base and completely remove it ③.

5. Store complete harness system for future use.


6. To remove the crotch buckle, please see section "Removing the crotch buckle".

7. Pull the up-front harness adjustment strap as far as it will go. This will lock the splitter plate into the base of the seat ④.

8. Roll strap tightly, secure roll with an elastic band and store rolled strap under the fabric seat pad in the opening of the up-front adjuster.

To reinstall, reverse the above steps.

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


To remove the crotch buckle:

1. Remove fabric crotch pad.
2. Push the buckle down through the fabric seat pad towards the plastic base of the child restraint.
3. Lift the fabric seat pad on the base of the child restraint to expose the buckle and slot ①.
4. Turn the metal plate so that it can slide through the opening in the plastic base ②. Continue to thread the crotch strap through the plastic base and completely remove it ③.

The buckle **MUST** be looped back through the seat base as shown ①.

! ENSURE that the strap is **COMPLETELY** in the slot and the crotch strap is **NOT** twisted.

⚠ WARNING!

FAILURE TO REINSTALL THE CROTCH STRAP
CORRECTLY COULD RESULT IN SERIOUS INJURY OR
DEATH.

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
5. Replace fabric seat pad and fabric crotch pad.

! ENSURE that the buckle is facing **FORWARD** ④.

PREPARING BOOSTER SEAT FOR USE

⚠ WARNING!

INCORRECT USE OF THE VEHICLE'S 3-POINT LAP-AND-SHOULDER SEAT BELT SYSTEM IS THE MOST COMMON MISTAKE PEOPLE MAKE. IT IS VERY IMPORTANT TO READ, UNDERSTAND AND FOLLOW ALL INSTRUCTIONS AND WARNINGS. DO NOT MISTAKE COMFORT FOR SAFETY. THE VEHICLE'S 3-POINT LAP-AND-SHOULDER SEAT BELT MUST BE SNUGLY ADJUSTED AT ALL TIMES.

! ENSURE the vehicle shoulder belt's position is correct for your child's size. **VERIFY OFTEN.**

! WARNING: DO NOT use this booster seat if the midpoint of your child's head is above the top of the booster seat headrest at its maximum height.

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⚠ WARNING!

FAILURE TO PROPERLY SECURE YOUR CHILD IN THE BOOSTER SEAT INCREASES THE RISK OF SERIOUS INJURY OR DEATH IN THE EVENT OF A SHARP TURN, SUDDEN STOP OR CRASH. In a crash, the booster seat provides **NO** protection for your child unless the vehicle 3-point lap-and-shoulder belt is properly secured. Follow these instructions carefully to properly secure your child.

! **WARNING:** You **MUST** use the vehicle 3-point lap-and-shoulder belt system when restraining the child in the booster seat ①.

! **WARNING:** Vehicle belts **MUST** be tight, adjusted to the child's body, and should not be twisted ①.

! **NEVER** use a vehicle lap-only belt to secure child in booster seat ②.

! **NEVER USE A LOOSE SHOULDER BELT.**

! **NEVER** position the shoulder belt under the child's arm ③ or place the shoulder belt behind the child's back ④.

! **Lap portion of lap-and-shoulder belt MUST be low and snug across the hips, with the pelvis firmly engaged, just touching the child's thighs.** The lap portion **MUST NOT** be lying across the child's stomach as this can cause major injuries to the child in case of a sharp turn, sudden stop or crash.

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! DO NOT allow the child to slouch, slide down in the seat, or create slack in the belt system **5**. If child will not keep vehicle belt properly positioned, **DO NOT** use this booster seat.

! DO NOT position vehicle belt over top of the armrests **6**.

! The vehicle seat headrest should not create a gap between vehicle seat and booster seat 7. If so, remove the vehicle headrest.

! DO NOT use any load bearing contact points other than those described in the instructions and marked in the child restraint.

IF IN DOUBT ABOUT THESE POINTS OR ANY QUESTIONS RELATED TO INSTALLATION, PLEASE CONTACT THOMASHILFEN.

INSTALLATION - BOOSTER SEAT

WARNING!

DUE TO VARYING ANGLES OF VEHICLE SEAT BELT RETRACTORS, EXCESS SLACK IN THE VEHICLE BELT SYSTEM MIGHT NOT BE REMOVED AUTOMATICALLY. Once the lap and shoulder portions of the seat belt are properly positioned and the buckle is attached, **FEED ANY REMAINING SLACK OF THE SEAT BELT MANUALLY BACK INTO THE RETRACTOR TO TIGHTEN THE WHOLE SYSTEM.**

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
1. Place child in the booster seat with child's back flat against the back of the booster seat.

2. Thread shoulder belt portion of the 3-point lap-and-shoulder belt system through the shoulder belt guide

① and ensure that shoulder belt lies across child's shoulder and chest, and **NOT** across child's face or neck.

3. Thread both sides of the lap belt and the lower end of the shoulder belt of the 3-point lap-and-shoulder belt system **UNDER BOTH** booster seat armrests ② and ensure that the belt is snug across child's hips. It is important to ensure that any lap strap is worn low down, so that the pelvis is firmly engaged, just touching their thighs and is **NOT** across their stomach.

4. Buckle the vehicle seat belt into place and ensure it is secure by gently pulling on it.

! BOTH SIDES OF THE LAP BELT AND THE LOWER END OF THE SHOULDER BELT MUST PASS UNDER BOTH ARMRESTS.

! Any straps holding the restraint to the vehicle should be tight, any straps restraining the child should be adjusted to the child's body, and all straps should not be twisted.

! Teach your child that seat belts are not toys and should **NEVER** be played with, regardless of whether the car is in motion or not.

! **ALWAYS** be aware of what your children are doing in the back seat.

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
TO INSTALL:

1. Slide cup holder into place over cup holder clip **1 2** on either side until locked in place.

TO REMOVE:

1. Pull cup holder up and outward.

**AIRCRAFT INSTALLATION**

This child restraint is **NOT CERTIFIED** for aircraft use. Aircraft seats are not equipped with 3-point lap-and-shoulder belts which are required for use with this child restraint.

CARE AND MAINTENANCE**CLEANING THE PLASTIC SHELL:**

- **SPONGE CLEAN** using warm water and mild soap.
- **TOWEL DRY.**
- **DO NOT** use solvents, abrasive cleaners, or disassemble child restraint in any way, except as instructed.

CLEANING HARNESS STRAPS:

- Spot clean with a damp cloth. **DO NOT BLEACH.** **DO NOT** immerse the straps in water. Doing so may weaken the straps.

CLEANING METAL AND PLASTIC HARNESS PARTS:

- Use mild soap and cool water on a damp cloth.

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GB

! NEVER lubricate the buckle. If buckle or harness adjuster sticks or will not remain tight, check for foreign objects such as food, dirt, etc. Remove objects with tweezers and/or clean with a cotton swab. If you cannot get the harness to remain tight with the buckle securely locked, **DO NOT USE** this child restraint. Call Thomashilfen to purchase a replacement harness or buckle.

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**INSTRUCTION MANUAL
FOR SPECIAL NEEDS
CAR SEAT**

Assembly and user instructions for special needs
version (supplement)

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- Gruppe I–III, 9-36 kg

Warnings	37
Intended Purpose / Indication.....	38
Application risks and contraindications.....	38
Symbols and warning signs	39
Delivery contents / Standard equipment	40
5-Point-Positioning belt	40
Head pads.....	41
Belt conector	42
Abduction block – pad (accessory).....	42
Seat wedge (accessory)	43
Seat wedge, inside (accessory).....	44
Table (accessory)	45
Lateral seat pads (accessory).....	46
Lateral side supports (accessory)	46
Crotch pad, large (accessory)	47
Adapter plate (accessory).....	47
Swivel plate (accessory)	49
Cushion for swivel base (accessory)	50
Back stabilizer (optimal extra)	51
Inclination plate (accessory).....	52
Footrest (accessory) and footrest, low (accessory)	53
Stand (accessory)	55
Seatfix adapter (accessory)	56
Stabilising bar (accessory)	57
Seat depth extension + 13 cm (accessory)	58
Abduction block, adjustable (accessory)	59
Seatfix-connector	60
Cleaning, disinfection and care instructions.....	61
Recycling/re-use	61
Maintenance schedule – annual inspection /re-use ..	61
Preparation	61

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Please read this and the assembly and user instructions provided with the seat first to ensure safe, easy and comfortable handling of your new car seat for you and especially your child.

These instructions are a supplement to the original assembly and user instructions. Some aspects in handling the car seat change due to modifying the car seat for use as special needs car seat. These changes are described in these supplementary instructions.

In this context, we would like to make you aware of our website www.thomashilfen.de where you will find up-to-date information and important product safety information or information on any product recalls.

If you still have any questions or experience any problems, please contact your competent retailer or us directly.

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Our address for customers in the USA:

ExoMotion LLC, dba Thomashilfen
7936 Occidental Ave S
Seattle WA 98108
USA
phone: 866 870 2122 (toll free in U.S.)
www.thomashilfen.us

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GENERAL INFORMATION

**WARNING!**

Based on the ECE R 44/04 requirements, the Major special needs car seat has successfully passed a stress test carried out by Dekra (Klettwitz). ECE Group I-III, 9 – 36 kg:

During this test it has been confirmed that the 5-point positioning belt, which is NOT safety-relevant, is suitable for use as additional positioning belt for children of ECE Group II-III (15kg-36kg). ALWAYS secure your child and the seat with the vehicle's 3-point safety belt and refer to the information regarding correct installation in the first part of the user instructions.

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The Major special needs car seat has been designed for the transport in vehicles of children with disabilities including the following conditions:

- infantile cerebral palsy
- hypotension of unknown origin
- muscular dystrophy
- various other syndromes

This car seat is characterised by its shell-like design and the provided pad system, which provides additional seating safety, stabilisation of head and pelvis, and support/control of the feet.

The Major special needs car seat therefore adjusts perfectly to the individual needs of children with disabilities.

APPLICATION RISKS AND CONTRAINDICATIONS

Attention: Please ensure that the clearance between seat and dashboard or front seat is at least 600 mm when using the Major special needs car seat in Group I.

Attention: Please ensure that the clearance to the roof is sufficient when using the Major special needs car seat for bigger children. This applies especially when using the adapter plate, inclination plate and/or the swivel plate.

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Specification labels / loading capacity stickers

The specification label of the car seat is attached at the bottom of the seat.

The specification label includes the name of the manufacturer, including the address; the type designation, the maximum loading capacity and the CE label.

**Bar code label**

The bar code label is attached under the seat.

The bar code label includes the item number, designation, serial number and production date of the product.

TECHNICAL DATA AND VERSIONS

Seat depth	32 cm
Seat width, front (excl./incl. pads)	31,5 – 21,5 cm
Seat width, back (excl./incl. pads)	26 – 19 cm
Back height, seat	62-74 cm
Shoulder height, strap guide	25,5 - 42 cm
Overall dimensions	56 x 71 x 40 cm
Weight of seat	approx. 7 kg
Max. user weight	9-36 kg (Group 1-3)
Isofix provided	yes

ITEM NUMBER / DESIGNATION

9800006100000	(610)	Major, grey
9800006110000	(611)	Major, red
9800006120000	(612)	Major, blue
9800006130000	(613)	Major, green
9800006140000	(614)	Major, magenta

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9800006100000	(610)	Major, grau
9800006110000	(611)	Major, rot
9800006120000	(612)	Major, blau
9800006130000	(613)	Major, grün
9800006140000	(614)	Major, magenta

The Major special needs car seat is delivered fully assembled and includes the following standard equipment:

- headrest
- crotch pad
- seat insert
- Seatfix connector

5-POINT POSITIONING BELT

General information on using the belt can be found in the first part of the instruction manual. As described in the general information, the 5-point positioning belt is NOT safety-relevant only in ECE Groups II-III.

Belt routing 3-point-belt

Important: Make sure that the 3-point seat belt is positioned in accordance with the front part of the instruction manual and the warning on the seat.

- 1 Warning on cushion

Attention: Make sure that the 3-point seat belt is positioned in accordance with the front part of the instruction manual and the warning on the seat.

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- 1 Warnhinweis am Polster

Pull the covers of the pads over the laterals of the headrest as tightly as possible. In addition, fix the head pads in place by using the elastic bands in the headrest. Push buttons in the seat shell provide additional stability.

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2 a: Druckknopf Kopfpelotte

2 b: Kopfpelotte befestigt

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To use the chest clip, unthread the 5-point positioning belt of the seat. The chest clip is then threaded through the ends of the 5-point positioning belt. Please make sure that the belts are not twisted and have been threaded properly. To finish, thread in the 5-point positioning belt again (in reverse order to pulling out).

③ : Chest clip in use

Attention: Do not use the chest clip for Group 1 (9-18 kg).

Attention: Information on correct unthreading and threading can be found in the first part of the instruction manual.

ABDUCTION BLOCK – PADDED (ACCESSORY) (ITEM NO.: 670)

The abduction block is fitted to the seat pad by using the belt strap on the seat insert. For this, pull the belt end through the openings in the pad and fasten again with the push button on the abduction block.

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SEAT WEDGE, BELOW (ACCESSORY) (ITEM NO.: 671)

The seat wedge is positioned centrally beneath the car seat. The high edge is at the front to achieve the ideal rest position for the child.

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

5 b: Endposition

The seat wedge (inside) can be placed on the seat in three different positions:

1. You achieve an active sitting posture by placing the seat wedge on the seat with the high edge facing backward and using the hook-and-loop fastener to fix in place.
2. For a more passive sitting position/rest position, attach the seat wedge with the flat side facing the back.
3. To reduce the depth of the seat or for a larger seat angle, attach the seat wedge upright at the back of the seat with the wide side facing downward.

Attach the two outer pads for additional protection at the sides.

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Place the table on both armrests on the left and right side. The table is fixed on both sides with two hook-and-loop fasteners each through the relevant D-rings in the seat cushion. In addition, the lower arm support must be fixed with a buckle through the sides of the back cushion.

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7 a: D-Ringe

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HIP SUPPORTS (ACCESSORY) (ITEM NO.: 674)

To attach the hip supports, simply pull the hoods over the armrests and, in addition, fix under the seat liner using the hook-and-loop fasteners. The seat insert can also be removed via the belt buckle and the openings in the cushion.

8 a: Hip supports over armrests

8 b: Attach under seat insert

LATERAL SIDE SUPPORTS (ACCESSORY) (ITEM NO.: 675)

Using hook-and-loop fasteners, the lateral side supports can be attached in various heights to the left and right in the back area. To attach properly, push the pads against the loop strap provided at the back of the seat.

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CROTCH PAD, LARGE (ACCESSORY) (ITEM NO.: 676)

The large crotch pad is pulled over the belt shield with the belt open.

Attention: The crotch pad from the standard equipment must be removed in this case.

10: Crotch pad over belt buckle

Attention: The crotch pad, large, cannot be used for Group 1 (9-18 kg).

ADAPTER PLATE (ACCESSORY) (ITEM NO.: 500)

To use the car seat with the swivel plate or the accessories below, mount the adapter plate to the seat first. First the supplied spacers must be attached to the seat. Insert the screw from above into the pre-drilled hole and then tighten the spacer from below. To do this, place the adapter plate with the tubes facing upwards on the base of the seat and fasten the plate with the screws provided in the cavity anchors and the mounted spacers.

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11 e: Tightening the screw

11 f: Mounting the spacer

Important: When using the adapter plate, the supplied spacers must always be mounted.

Important: If additional accessories are used, please make sure that they can only be pulled out until the marking is visible.

Attention: The modular construction of the plate system must always be carried out in the following order (from the top):

1. Car seat
2. Adapter plate
3. Tilting base
4. Swivel base

SWIVEL BASE (ACCESSORY) (ITEM NO.: 501)

The swivel base is supplied fully assembled and must be attached to the adapter or tilting base using the hook-and-loop fasteners at the top. The direction of rotation can be changed if necessary. For this purpose, the locking bolt, the metal web and the rotary screw must be mounted on the other side. In addition, the turn-limiter must be turned by 180°.

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- 12 a: Removal of the locking bolt
- 12 b: Removal of the metal web
- 12 c: Turn-limiter clockwise
- 12 d: Turn-limiter counter clockwise

Please note: We recommend that the plates be placed from back to front.

Please note: To enable a larger angle of rotation, mount the locking bolt on the other side without changing the turn-limiter.

Attention: The swivel base cannot be used for Group 1 (9-18 kg).

Attention: The swivel base can only be used in connection with an adapter plate.

Attention: In a modular construction, the swivel base must always be at the bottom.

CUSHION FOR SWIVEL BASE (ACCESSORY) (ITEM NO.: 513)

The cushion is placed centrally under the swivel base, with the wider side in front so that it matches the contour of the swivel base. It ensures that the tubes of the swivel base do not lie directly on the vehicle seat and the child car seat is brought into a higher sitting position, which might improve the swivel function.

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BACK STABILIZER (ACCESSORY) (ITEM NO.: 66308)

Please note: We recommend to mount a back stabilizer to the Major when a turning plate is used to prevent the backrest from accidentally falling down during turning.

The back stabilizer is to be put at the seat and needs to be screwed into the cavity anchors..

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TILTING BASE (ACCESSORY) (ITEM NO.: 502)

The height adjustment plates must each be attached by fixing the two supplied screws to the bottom plate. The tilting base must then be attached to the adapter plate by using the hook-and-loop fasteners at the top.

15 a: Hook-and-loop fasteners on the tilting base

15 b: Tilting base arrangement

Please note: We recommend to put on the plate from back to the front.

Attention: The tilting base cannot be used for Group 1 (9-18 kg).

Attention: The tilting base can only be used in connection with an adapter plate.

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**FOOTREST, SIZE 1 (ACCESSORY)
(ITEM NO.: 503)
AND FOOTREST, SIZE 2 (ACCESSORY)
(ITEM NO.: 504)**

First, loosen the locking screw at the front of the adapter plate's outer tube mounting and insert the guide tubes of the footrest as far as possible into the mounting provided beneath the adapter plate. Now tighten the locking screw beneath the adapter plate again.

To stabilise the stand on the vehicle floor, adjust the height of the supporting rods by screwing them in or out until the rods touch the floor of your vehicle. Now lock the supporting rods with the star knobs to prevent unintentional unscrewing. For this, turn the star knob screws upward as far as possible and tighten securely.

Attention: The adapter plate or child car seat should always be just about on the vehicle's seat!

Lower leg length:

The lower leg length is adjusted by changing the height of the footrest. For this, loosen the two outer Allen screws and tighten them again after adjustment.

Foot angle:

Loosen the four screws left and right at the bottom of the footrest to adjust the desired foot angle with the help of the guides. Then tighten the screws again.

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16 c: height adjustment

16 d: Securing of the foot rest

Attention: If the tilting and/or the swivel base is used, then the footrests must always be set at the lowest position in the tubes.

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Note: Using the stand can improve the stability of the entire car seat. The stand can be used when the child car seat is stabilised by the Seatfix adapter and a footrest is not provided.

First, loosen the locking screw at the front of the adapter plate's outer tube mounting and insert the guide tubes of the stand as far as possible into the mounting provided beneath the adapter plate. Now tighten the locking screw beneath the adapter plate again.

To stabilise the stand on the vehicle floor, adjust the height of the supporting rods by screwing them in or out until the rods touch the floor of your vehicle. Now lock the supporting rods with the star knobs to prevent unintentional unscrewing. For this, turn the star knob screws upward as far as possible and tighten securely.

17: Securing of the stand

Attention: The adapter plate or child car seat should always be just about on the vehicle's seat!

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Note: Requirement for using the Seatfix adapter is that ISOFIX equipment is provided in the vehicle (mostly at the back on the passenger side). This must be checked prior to purchase!

First, loosen the locking screws at the back of the adapter plate's outer tube mounting. Then place the Seatfix adapter on the vehicle's seat. The hole at each end of the guide tubes must be facing downward to the vehicle's seat. Insert the Seatfix connectors (with yellow release buttons) in the Isofix anchors of the vehicle. Pull and check whether the connectors have engaged properly. Now insert the Seatfix adapter in the mounting at the back of the swivel plate (where the plate is rounded) with the locking screw loosened.

Adjust the extension length of the Seatfix adapter to your vehicle. When attached to the swivel plate, the child car seat must be in contact with the backrest of the vehicle's seat. Fix the Seatfix adapter now with the loosened locking screw.

Please note: When using the Seatfix adapter, additional stability can be provided by supporting the swivel plate / seat at the front toward the floor of the vehicle. For this, use the footrest (Item no. 503 or 504) or the stand (Item no. 505).

Attention: The Seatfix adapter for the swivel plate does not replace the safety belt, but is used for positioning the car seat. The child / seat must still be secured with the 3-point belt of the vehicle in accordance with the instruction manual.

Attention: If the inclination and/or the swivel plate is used, then the Seatfix adapter must always be set at the lowest position in the tubes.

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18 b: built-in seatfix adapter

STABILIZING BOW (ACCESSORY) (ITEM NO.: 507)

First, loosen the locking screw at the back of the tube mounting beneath the adapter plate. Insert the stabilizing bar into the tube mountings at the rounded part at the back.

Now place the adapter plate on the seat so that the stabilising bar is inserted as far as possible between the back and the seat of the vehicle.

Following this, adjust the extension length of the stabilising bar to your vehicle. When attached to the swivel plate, the child car seat must be in contact with the backrest of the vehicle's seat.

Now tighten the locking screw properly.

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Note: The stabilising bar has two sides which can be used. Try in your vehicle which side provides the most stability for the swivel plate.

Attention: If the inclination and/or the swivel plate is used, then the stabilising bar must always be set at the lowest position in the tubes.

SEAT DEPTH EXTENSION + 13 CM (ACCESSORY) (ITEM NO.: 508)

Using the guide tubes, the seat depth extension for the base is inserted as far as possible into tubes no. 2 and 4 and then fixed in position by tightening the locking screws beneath the plate.

20: Guide tube for seat depth extension

Please note: If you use a footrest (art.no.:503/504) or the stand (art.no.:505) we recommend to mount the seat depth extension first.

Please note: We recommend to route the adjustment strap below the seat depth extension.

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ABDUCTION BLOCK, ADJUSTABLE (ACCESSORY) (ITEM NO.: 509)

The abduction block for the base is inserted with the guide tube into the middle mounting and is tightened with the locking screw beneath the plate. To adjust the position of the abduction block to the specific requirements and size of the child, the depth is adjustable. Loosen the locking screw in the mounting and adjust the abduction block accordingly.

21 a: Abduction block guide tube

21 b: Abduction block locking screw

Attention: The adjustable abduction block cannot be used for Group 1 (9-18 kg).

Please note: If you use a footrest (art.no.: 503/504), the stand (art.no.: 505) or the seat depth extension (art.no.: 508) we recommend to mount the adjustable abduction block first.

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SEATFIX - CONNECTOR (ITEM NO.: 66289/66303)

Note: Requirement for using the Seatfix adapter is that ISOFIX equipment is provided in the vehicle (mostly at the back on the passenger side). This must be checked prior to purchase!

The Seatfix connectors are located beneath the seat shell of your child car seat. Prior to placing the child car seat in the vehicle, expand the connectors so they become visible. After placing the child car seat in the vehicle, push the connectors into the Isofix brackets until a click can be heard and both connectors display the green mark. Pull and check whether the connectors have engaged properly.

22 a: Expanding the Seatfix connectors

22 b: Expanded connectors

Attention: The Seatfix connectors can no longer be used and the Seatfix adapter (Item no. 506) must be used instead as soon as the adapter plate is used.

Attention: The Seatfix connector does not replace the safety belt, but is used to position the car seat. The child / seat must still be secured with the 3-point positioning belt of the vehicle in accordance with the instruction manuals.

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PADS

All pads are easy to remove. The 100% Trevira CS fabric is of high quality and washable at temperatures of up to 40°C. Wash pads separately and cover the hook-and-loop fasteners prior to washing. Please note that even high-quality fabrics can fade due to constant and intense exposure to sunlight or frequent washing. The pads are flame retardant (EN 1021-2).

HOOK-AND-LOOP FASTENERS

To maintain the functionality of the hook-and-loop fasteners, clean occasionally using a brush. When washing, always close or cover the hook-and-loop fasteners.

RECYCLING/ RE-USE

If the seat is no longer usable due to its condition, please take it to a competent local waste disposal company for recycling.

MAINTENANCE SCHEDULE –
ANNUAL INSPECTION / RE-USE

PREPARATION

An instruction manual must be provided with the product. If not available, please contact the manufacturer and request an instruction manual. Familiarise yourself with the functions of the product. If you don't know the product, read the instruction manual prior to the inspection. Clean the product prior to inspection.

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Member of the International
Support Association for the
Rehabilitation of Children and Youths e.V.

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Thomas Hilfen für Körperbehinderte GmbH & Co. Medico KG · Walkmühlenstr. 1 · D - 27432 Bremervörde · Germany